

**SEAT ROTATES  
360° AND LOCKS  
EVERY 90°!**



**Bath Stool with Padded Rotating Seat**

RTL12061 4 retail boxes/cs

**FEATURES:**

- The rotating padded seat reduces twisting and reaching while in the shower.
- The rotating top of the seat can be removed for easy cleaning.
- Comes with removable tray for storing personal items.
- Ideal for individuals who have balance or mobility issues.
- Seat height can be adjusted from 16.14 to 22.8.
- Weight capacity: 300 lbs.
- Limited Lifetime Warranty.



**RETAIL  
PACKAGED!**

**K.D. Teak Adjustable  
Height Bath Stool**

12350KDR-1 Triangular Seat, Knock Down, Tool Free Assembly, Retail, (Will become RTL12350KDR in 2011), 1/cs

12351KDR-1 Rectangular Seat, Knock Down, Tool Free Assembly, Retail, (Will become RTL12351KDR in 2011), 1/cs

10389G-1 Replacement Tips for 12350, 12351, 4/bx

**FEATURES:**

- Attractive high-quality teak is renowned for its durability and weathering properties.
- **Ideal for use in shower stall.**
- Teak will naturally take on a handsome silver grey patina after several months of exposure to water.
- Ideal for a woman's grooming stool.
- Lightweight, durable, adjustable legs.
- Large rubber tips reduce slipping.
- Limited Lifetime Warranty.



12350KDR-1



12351KDR-1



**TOOL FREE  
ASSEMBLY**



**STAINLESS STEEL THUMB  
SCREWS WITH PLASTIC WING  
NUTS ASSEMBLE TOOL FREE  
INTO STAINLESS STEEL INSERTS  
UNDER BLOW MOLDED SEAT.**



**SPECIFICATIONS**

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
RTL12061	15"	13.75"	15"	13.75"	22.8"	16.14"	7.4 lbs.	33 lbs.	300 lbs.
12350KDR-1	22"	22.5"	22"	16"	19"	13"	9 lbs.	11 lbs.	300 lbs.
12351KDR-1	17"	20"	17"	10"	18"	13"	7.5 lbs.	9 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

bath safety

To order call toll free: 877.224.0946

drive