



ARM WIDTH IS ADJUSTABLE!

RTL12087

Toilet Safety Rail

RTL12087 2 retail boxes/cs



FEATURES:

- Allows individuals to safely sit or stand.
- Easy assembly and installation.
- Comes with soft, padded handles for comfort.
- Handles are width adjustable from 16" to 18".
- Limited Lifetime Warranty.



COMPARE TO CAREX® #B368



RTL12079

Stand Alone Toilet Safety Rail

RTL12079 2 retail boxes/cs



FEATURES:

- Free standing around toilet.
- **Easy, tool free assembly using wing nuts.**
- Fits standard and elongated toilets.
- Comes with padded armrests for added comfort.
- Comes with a magazine rack for added convenience.
- Limited Lifetime Warranty.



TOOL FREE ASSEMBLY!

SPECIFICATIONS

Item #	Width Between Arms (smallest)	Width Between Arms (largest)	Depth	Height Seat	Product Weight	Carton Shipping Weight	Weight Limit
RTL12087	19.5"	22.5"	16.25"	11.5"	10.5 lbs.	24 lbs.	300 lbs.
RTL12079	NA	26.3" (overall width)	19.5"	NA	9.0 lbs.	25.3 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.