



**EURO-STYLE  
CLIP FOR  
HEIGHT  
ADJUSTMENT!**



12350KDR-1



**TOOL FREE  
ASSEMBLY**



12351KDR-1



**STAINLESS STEEL THUMB  
SCREWS WITH PLASTIC WING  
NUTS ASSEMBLE TOOL FREE  
INTO STAINLESS STEEL INSERTS  
UNDER BLOW MOLDED SEAT.**

### K.D. Euro-Style Bath Stool

- 12057\* Without Hygiene Cutout, 1/cs
- 12059\* With Two Hygiene Cutouts, 1/cs
- 12021C Replacement Euro-Style Clip, 1/ea

#### FEATURES:

- Two models available with hygiene cutouts.
- Easy to install legs without tools.
- Euro-style clip locks leg to seat
- Extra large suction cups provide added safety.
- Drainage holes prevent water from pooling on seat.
- Manufactured with 1" aluminum tubing to ensure lightweight stability.
- Impact resistant, composite seat is crack-proof and tarnish resistant.
- Limited Lifetime Warranty.

### K.D. Teak Adjustable Height Bath Stool

- 12350KDR-1 Triangular Seat, Knock Down, Tool Free Assembly, Retail, 1/cs
- 12351KDR-1 Rectangular Seat, Knock Down, Tool Free Assembly, Retail, 1/cs
- 10389G-1 Replacement Tips for 12350, 12351, 4/bx

#### FEATURES:

- Attractive high-quality teak is renowned for its durability and weathering properties.
- **Ideal for use in shower stall.**
- Teak will naturally take on a handsome silver grey patina after several months of exposure to water.
- Ideal for a woman's grooming stool.
- Lightweight, durable, adjustable legs.
- Large rubber tips reduce slipping.
- Limited Lifetime Warranty.



\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

#### SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12057/9	15"	15"	15"	15"	21.6"	17.7"	3.5 lbs.	5 lbs.	400 lbs.
12350KDR-1	22"	22.5"	22"	16"	19"	13"	9 lbs.	11 lbs.	300 lbs.
12351KDR-1	17"	20"	17"	10"	18"	13"	7.5 lbs.	9 lbs.	300 lbs.

bath safety  
To order call toll free: 877.224.0946

drive  
WWW.DRIVEMEDICAL.COM