

Wallaby Pediatric Folding Wheelchair

1P = 1-Pack Available

***WC19 Bus Transit Approved**

Model #:
Sold 1 each

1P **WB 1200**
1P **WB 1400**

Wallaby Pediatric Folding Wheelchair, 12", Red
Wallaby Pediatric Folding Wheelchair, 14", Blue



WB 1200



WB 1400

Accessories

WB 8021 Headrest Extension
Black nylon, width adjustable headrest extension provides additional head support for taller users. Will prevent head whiplash during rear end collision on a bus.

WB 8030 "H" Harness with Padded Reversible Strap Covers
Five-point harness with padded reversible strap covers increase users comfort and reduce friction. It features imitation lamb's wool on one side and velour on the other.

Weight Capacity 150 lbs.
80 lbs. when used in bus transport

Warranty Limited Lifetime

*The Wallaby is crash tested and meets bus transit requirements when used together with the headrest extension and "H" harness.

Features

- Powder coated colored steel frame with black cross brace.
- Easy to clean lightweight durable nylon upholstery.
- Height adjustable push handles.
- Flip-back padded desk length arms.
- Pelvic belt with reversible cover.
- 22" rear and 6" front swivel solid rubber tires mounted on composite wheels.
- Push-to-lock brakes
- Swing-away footrests with plastic plates and calf straps.
- Limited Lifetime Warranty

Specifications

| Item # | Overall Width (open) | Overall Width (closed) | Seat Depth | Seat Width | Seat to Floor (front) | Seat to Floor (rear) | Seat to Footrest | Back Cane Height | Push Handle Height |
|---------|----------------------|------------------------|------------|------------|-----------------------|----------------------|------------------|------------------|--------------------|
| WB 1200 | 20" | 12" | 12" | 12" | 19" | 18" | 13" - 16" | 14.5" | 33" - 39" |
| WB 1400 | 22" | 12" | 12" | 14" | 19" | 18" | 13" - 16" | 15" | 33" - 39" |

| Item # | Armrest Height | Overall Length w/Footrests | Overall Length w/o Footrests | Back Angle | Product Weight w/Footrests | Product Weight w/o Footrests | Weight Capacity | Weight Capacity When Used in Bus Transport |
|---------|----------------|----------------------------|------------------------------|------------|----------------------------|------------------------------|-----------------|--|
| WB 1200 | 7" | 37" | 28.5" | 90° | 37 lbs. | 33 lbs. | 150 lbs. | 80 lbs. |
| WB 1400 | 7" | 37" | 28.5" | 90° | 37 lbs. | 33 lbs. | 150 lbs. | 80 lbs. |



WB 1400, WB 8021, WB 8030

To order call toll free: 877.224.0946

WENZELITE

Wenzelite ReRehab **drive**