



**New Neo-Snap Connector!**



**Limited-reuse BP cuff**

**GE CRITIKON NEONATAL**

**Blood Pressure Cuffs**

Providing big benefits for your tiniest patients, GE Healthcare's neonatal cuffs are designed with the accuracy you trust and the gentle comfort you expect. We offer a full range of neonatal sizes with the choice of extra soft, absorbent material or easy-to-clean vinyl. The playful, color-coded animal designs help you find the right size at a glance and they are economically priced for single-patient use.



**A trusted choice for clinical excellence**

- Tested for performance accuracy and clinically validated to help you deliver quality patient care
- Full range of neonatal sizes in playful, color-coded animal designs to help you find the right size at a glance, enabling more accurate measurement
- Choice of extra soft, absorbent material (SOFT-CUF\*) or easy-to-clean vinyl (CLASSIC-CUF\*)

**Constructed for patient safety**

- Limited reuse, yet economical enough for pairing with a single baby throughout their hospital stay to help with your infection control efforts
- Cuff materials are not made with natural rubber latex or DEHP, making them safer for the babies
- New Neo-Snap\* connector meets international standards, prevents misconnections to luer fittings as well as helps prevent accidental disconnections for patient safety

**Designed to improve workflow**

- New Neo-Snap\* connector does not have to be aligned to make a secure connection, enabling you to focus more on the baby
- Cuffs compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- Full range of adapters and connectors are compatible with automated monitors from GE, Spacelabs\*\*, Philips\*\* and Datascope\*\*, plus on a variety of manual devices, including wall manometers, enabling a universal connection system throughout the facility

**Cuff assortment packs available**

- For low-volume cuff utilizers, several types of assortment packs available for easier ordering and more optimal inventory management



Neonatal			GE Healthcare				Philips/HP**, Siemens** Spacelabs and Datascope			
										
			2-Tube Luer Connector		2-Tube Neo-Snap Connector		1-Tube Luer Connector		1-Tube Neo-Snap Connector	
Size	Color	Limb Circumference	SOFT-CUF	CLASSIC-CUF	SOFT-CUF	CLASSIC-CUF	SOFT-CUF	CLASSIC-CUF	SOFT-CUF	CLASSIC-CUF
1	Orange	3-6 cm	2521	2638	SFT-N1-2B	CLA-N1-2B	2121	2138	SFT-N1-1B	CLA-N1-1B
2	Lt. Blue	4-8 cm	2422	2633	SFT-N2-2B	CLA-N2-2B	2122	2133	SFT-N2-1B	CLA-N2-1B
3	Green	6-11 cm	2523	2628	SFT-N3-2B	CLA-N3-2B	2124	2128	SFT-N3-1B	CLA-N3-1B
4	Navy	7-13 cm	2524	2623	SFT-N4-2B	CLA-N4-2B	2125	2123	SFT-N4-1B	CLA-N4-1B
5	Rose	8-15 cm	2525	2619	SFT-N5-2B	CLA-N5-2B	2126	2119	SFT-N5-1B	CLA-N5-1B

Neonatal Assortment Pack	Cuff Type	Part Number	Description
2-Tube Luer Connectors	SOFT-CUF	2694	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
	CLASSIC-CUF	2693	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
2-Tube Neo-Snap Connectors	SOFT-CUF	2059305-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
	CLASSIC-CUF	2059306-002	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
1-Tube Luer Connectors	SOFT-CUF	2696	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
	CLASSIC-CUF	2680	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
1-Tube Neo-Snap Connectors	SOFT-CUF	2059305-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)
	CLASSIC-CUF	2059306-001	2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5 (20 cuffs/pack)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [cuffs.gehealthcare.com](http://cuffs.gehealthcare.com)

GE Healthcare  
P.O. Box 900, FIN-00031 GE, Finland  
GE Direct United Kingdom: +44 (0)800 0329201

[www.gehealthcare.com](http://www.gehealthcare.com)



GE imagination at work

© 2013 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE, GE Monogram are trademarks of General Electric Company.

\*Trademarks of General Electric Company.

\*\*All third party trademarks are the property of their respective owner.

EMEA DOC1313640 04/13  
(Global version DOC1094824 rev 2)