

## Grab Bars

- Heavy duty steel construction with rust resistant finishes
- Available in chrome plated, brushed nickel, or white enamel finish
- Knurling provides a good gripping surface
- Mounts vertically, horizontally or diagonally 1 1/2" from wall
- Attractive retail packaging
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0241
- Lifetime limited warranty

### Brushed Nickel Finish:

|              |     |      |
|--------------|-----|------|
| <b>3012A</b> | 12" | 1 ea |
| <b>3016A</b> | 16" | 1 ea |
| <b>3018A</b> | 18" | 1 ea |
| <b>3019A</b> | 24" | 1 ea |

### Chrome Plated Finish:

|              |     |      |
|--------------|-----|------|
| <b>2012A</b> | 12" | 1 ea |
| <b>2016A</b> | 16" | 1 ea |
| <b>2018A</b> | 18" | 1 ea |
| <b>2019A</b> | 24" | 1 ea |
| <b>2021A</b> | 32" | 1 ea |

### White Versaguard Finish:

|               |     |      |
|---------------|-----|------|
| <b>40033A</b> | 12" | 3/cs |
| <b>40433A</b> | 16" | 3/cs |
| <b>40133A</b> | 18" | 3/cs |
| <b>40233A</b> | 24" | 3/cs |
| <b>40333A</b> | 32" | 3/cs |

