

Platinum Collection 3-in-1 Steel Padded Drop Arm Commode

- Arms drop independently below seat level for safe lateral transfers
- Easy-to-release drop arm mechanism may be operated by patients lacking manual dexterity
- Arms automatically lock in raised position
- Watertight, heat-sealed padded seat is easy to clean and comfortable
- May be used as a bedside commode, raised toilet seat or toilet safety frame
- Welded steel construction increases strength and stability
- Chrome-plated steel frame is scratch resistant and easy to clean
- Cushioned armrests enhance user comfort
- Seat height readily adjusts from 19 ½" to 23 ½", in 1 inch increments
- Comes complete with commode pail, cover and splash shield
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0165
- Limited lifetime warranty

6437A 1 ea



Specifications

Width Between Arms	18"
Overall Width	23"
Seat Width	14"
Seat Depth.....	16 ½"
Seat Height Range	19 ½" - 23 ½"
Maximum Weight Capacity.....	300 lb EVENLY DISTRIBUTED
Replacement Tip.....	9130



GRAHAM-FIELD

www.grahamfield.com