

Bariatric Bath Seat

- Durable blow-molded plastic seat with optional backrest
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- Comes assembled with tool-free back attachment
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0245
- Lifetime limited warranty

without backrest

7932A 2/cs

with backrest

7939A 2/cs

7939A-1 1 ea

Specifications	7932A	7939A
Overall Width.....	18"	18"
Overall Depth.....	16 ¾"	20 ½"
Seat Width.....	17 ¾"	17"
Seat Depth.....	12"	16 ¼"
Minimum Seat Height.....	14"	16"
Maximum Seat Height.....	17"	20"
Height of Back Support Above Seat.....	N/A.....	13"
Maximum Weight Capacity.....	600 lb EVENLY DISTRIBUTED	

