

SURELIFT™ Universal Slings with Full Head Support

- Provides excellent patient comfort and security
- Reinforced back for support
- Non-mesh slings have padding on legs and back
- Height-Adjustable straps
- Easy to use
- Machine washable
- Six point connectivity is designed to enhance positioning and safety
- Designed to meet the requirements of HCPCS Code: E0621
- 6 month limited warranty

		Body Width	Body Length	Total Width	Total Length	Weight Capacity
UPH862	Medium	38"	34"	78"	75"	400 lb
UPH863	Large	44 ½"	40"	80"	82"	400 lb
UPH864	XX-Large, Bariatric	45"	44"	85"	100"	600 lb

Also available in mesh:

UMH832	Medium	38"	34"	78"	75"	400 lb
UMH833	Large	44 ½"	40"	80"	82"	400 lb

