



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Three Drawer Mini Line Infection Control Cart (Short Cabinet), Key-locking, Standard Package #3143K

Specifications



Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 24 vertical inches (61 cm) drawer space
- Key-locking system with 2 keys
- 3" (7.6 cm) single-wheel stem caster set, full swivel – two with brakes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Pre-threaded accessory mounting holes
- Latex free
- Twelve year warranty
- Cart requires caster installation and accessory assembly after delivery
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) - 6-inch (15.2 cm)
- (2) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32 cm x 42.5cm)

Assembled Cart Dimensions

- 34.4" H x 18" W x 18"D
(87.4 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 71 pounds (32.2 kg)
- Approximate shipping weight: 90 pounds (41 kg)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black