



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Five Drawer Mini24 Anesthesia Cart with Bumper and 5" Wheels, Key Lock #4245K

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices

Specifications



Color Shown is Sand.
Additional Colors Available.

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 24 vertical inches (61 cm) of drawer space
- Key lock with two keys
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5" (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

- 18.5" W x 16.75" D (47 cm x 42.5cm)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Assembled Cart Dimensions:

- 36.5" H x 29.25" W x 22" D
(92.5 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 103 pounds (46.8 kg)

Shipping Dimensions:

- 49" H x 33" W x 26" D
(124.5 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 142 pounds (64.4 kg)



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Five Drawer Mini24 Anesthesia Cart with Bumper and 5" Wheels, with Basic Electronic Lock #4245E

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices

Specifications



Color Shown is Sand.
Additional Colors Available.

Shipping Dimensions:

- 49" H x 33" W x 26" D
(124.5 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 142 pounds (64.4 kg)

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 24 vertical inches (61 cm) drawer space
- Basic electronic pushbutton lock with keylock override
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5" (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

- 18.5" W x 16.75" D (47 cm x 42.5cm)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Assembled Cart Dimensions:

- 36.5" H x 29.25" W x 22" D
(92.5 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 103 pounds (46.8 kg)



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Five Drawer Mini24 Anesthesia Cart with Bumper and 5” Wheels, Electronic lock cart with keypad access #4245EKC

Harloff’s Mini24 Line offers a 24” wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices



Color Shown is Hammertone Blue
 Additional Colors Available.

Assembled Cart Dimensions:

- 36.5” H x 29.25” W x 22” D
 (92.5 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 106 pounds (48.2 kg)

Shipping Dimensions:

- 49” H x 33” W x 26” D
 (124.5 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 158 pounds (71.7 kg)

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 24 vertical inches (61 cm) of drawer space
- Electronic lock cart with keypad access and key lock override with 2 keys
- Drawers automatically re-lock after access (pre-set or programmable time)
- Access rights may be assigned for up to 250 users per cart
- Allows for supervisor assignments; dual credential access available
- May add, edit, or delete users and user groups
- Battery-operated with six D-cell batteries; life of one year (minimum) at 25 to 50 uses/day
- Lock has a LED low battery indicator
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5” (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

- 18.5” W x 16.75” D (47 cm x 42.5cm)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Five Drawer Mini24 Anesthesia Cart with Bumper and 5” Wheels, Electronic lock with proximity card and keypad access #4245ELP

Harloff’s Mini24 Line offers a 24” wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices



Color Shown is Hammertone Green
 Additional Colors Available.

Assembled Cart Dimensions:

- 36.5” H x 29.25” W x 22” D
 (92.5 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 106 pounds (48.2 kg)

Shipping Dimensions:

- 49” H x 33” W x 26” D
 (124.5 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 158 pounds (71.7 kg)

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 24 vertical inches (61 cm) of drawer space
- Electronic lock with proximity card and keypad access, and key lock override with 2 keys
- Drawers automatically re-lock after access (pre-set or programmable time)
- Access rights may be assigned for up to 250 users per cart
- Allows for supervisor assignments; dual credential access available
- May add, edit, or delete users and user groups
- Battery-operated with six D-cell batteries; life of one year (minimum) at 25 to 50 uses/day
- Lock has a LED low battery indicator
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5” (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

Internal Drawer Dimensions:

- 18.5” W x 16.75” D (47 cm x 42.5cm)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)