



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

## Five Drawer Mini24 Isolation Cart with Bumper and 5" Wheels with Basic Electronic Lock #4255E

Harloff's Mini24 Line offers a 24" wide cabinet with more internal storage space than our standard Mini Line. Select the cart best suited for your needs from our multiple cabinet heights, drawer configurations, and locking system choices.

### Specifications

#### Construction:

- Painted steel cabinet and drawers

#### Features:

- Space-saver design allows supply storage in areas with space limitations
- Cabinet has 30 vertical inches (76.2 cm) of drawer space
- Basic electronic pushbutton lock with keylock override
- Protective bumper extends wheel base and stabilizes cart when drawers are opened
- 5" (12.7 cm) Premium caster set, full swivel – two with brakes
- Ships with wheels and bumper assembled
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty on manufactured parts
- Durable hard-baked powder coat finish available in seventeen colors

#### Internal Drawer Dimensions:

- 18.5" W x 16.75" D (47 cm x 42.5cm)

#### Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (3) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

#### Assembled Cart Dimensions:

- 42.75" H x 29.25" W x 22" D (108.4 cm H x 74.1 cm W x 55.9 cm D)
- Cart Net Weight: 119 pounds (54.1kg)



Color Shown is Hammertone Blue  
 Additional Colors Available.

#### Shipping Dimensions:

- 57" H x 33" W x 26" D (144.8 cm H x 83.8 cm W x 66 cm D)
- Approximate shipping weight: 156 pounds (70.76 kg)