



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**3 Drawers with Storage Compartment, GP Line, Key Locking,
 Standard Package
 #6150**



Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Key lock with four keys
- 18" Storage compartment with locking door, and adjustable shelf
- Narrower cabinet for limited space
- Full wrap-around bumper on separate stabilizing frame
- 5" (12.7 cm) ball-bearing premium casters, all full swivel-one tracking-one braking
- Aluminum push handle / side rail
- Ball-bearing, full extension drawer slides
- Cart is delivered fully assembled
- Latex free
- Compatible with select accessories and drawer dividers
- Twelve year warranty
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (2) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 18" (45.7 cm) Storage compartment

Sixteen available color choices:

***Antimicrobial Finish Available**

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Internal Drawer Dimensions:

- 20" W x 17.25" D (50.8 cm W x 13.8 cm D)

Assembled Cart Dimensions

- 40.25" H x 28.5" W x 22.25" D (102.0 cm H x 72.4 cm D x 56.5 cm W)
- Cart weight: 125 pounds (56.8 kg)
- Shipping weight: 140 pounds (63.5 kg)