



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## 3 Drawers with Storage Compartment, GP Line, Breakaway Locking, Specialty Package #6161



### Specifications

#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Stainless steel break-away lock bar
- Pack of 50 break-away seals
- 18" Storage compartment with locking door, and adjustable shelf
- Narrower cabinet for limited space
- Full wrap-around bumper on separate stabilizing frame
- 5" (12.7 cm) ball-bearing premium casters, all full swivel-one tracking-one braking
- Stainless steel top rail / push handle
- Removable plastic top
- Ball-bearing, full extension drawer slides
- Cart is delivered fully assembled
- Latex free
- Compatible with select accessories and drawer dividers
- Twelve year warranty
- Durable, hard-baked powder coat finish

#### Drawer Configuration:

- (2) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 18" (45.7 cm) Storage Cabinet

#### Internal Drawer Dimensions:

- 20" W x 17.25" D (50.8 cm W x 13.8 cm D)

#### Sixteen available color choices:

##### \*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

#### Assembled Cart Dimensions

- 43" H x 28.5" W x 22.25" D  
(109.2 cm H x 72.4 cm D x 56.5 cm W)
- Cart weight: 125 pounds (56.9 kg)
- Shipping weight: 140 pounds (63.5 kg)