



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Long Term Care Treatment Cart – Standard Package #6220

Specifications

Construction:

- Painted steel

Features:

- Tall cabinet
- Dual flush-folding full height utility doors with slam lock feature
- Corner bumpers
- 5" (12.7cm) sealed ball-bearing premium casters, all full swivel – one with brake, one directional
- 20 standard size storage bins – 5.5" W x 9.5" D x 2.5" H (14 cm x 24.1 cm x 6.4 cm) located on 5 fixed shelves (4 bins per shelf)
- One roll-out shelf at bottom with deep storage tub
- Tub dims: 22.38" W x 17.38" D x 8" H (56.9 cm x 44.2 cm x 20.4 cm)
- Removable plastic top
- Cart is delivered fully assembled
- Certified "Latex Free"
- Hard-baked powder coat finish
- Five year warranty



Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
Light Blue, Navy, Light Gray,
Mauve, Red, Teal, Yellow,
Hammer Tone Blue, Hammer Tone Gray
and Hammer Tone Green

Assembled Cart Dimensions:

- 41.5" H x 32" W x 22" D
(105.4 cm x 81.3 cm x 55.9 cm)
- Cart Net Weight: 147 lbs (66.9 kg)
- Approx. Shipping Weight: 182 lbs (82.8 kg)