

Three Drawer Infection Control Cart with Key Lock Standard Package #6510

Specifications



Construction:

- All steel

Features include:

- Tall cabinet with key locking system
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) premium rubber single-wheel caster set, full swivel – one with brake, one directional
- Ball-bearing drawer slides
- Double-walled drawer faces with no rough edges
- Removable plastic top
- Aluminum push handle
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish

Twelve available color choices:

- Yellow (-YL), Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Drawer configuration:

- 1 x 6-inch (15.2 cm)
- 2 x 12-inch (30.5 cm)

Assembled cart dimensions:

- 41.5" H x 32" W x 22" D (105 cm H x 81.3 cm W x 55.9 cm D)
- Cart weight: 140 pounds (63.5 kg)
- Approx. shipping weight: 175 pounds (79.6 kg)