



650 Ford Street  
Colorado Springs, CO 80915  
Toll:(800)433-4064 \* Ph: (719)637-0300  
Fax:(719)597-8273 \* www.harloff.com

## 4 Drawer Workstation (Short) with Basic Electronic Pushbutton Lock, Standard Package #7350E

### Specifications

#### Construction:

- Sturdy, uni-body steel construction

#### Features:

- Basic Electronic Pushbutton Lock
- Key Lock Override
- Factory Set User and Supervisor Code
- Cabinet may be programmed in the field
- 24" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Pre-drilled threaded holes for side mounted accessories
- Full wrap-around bumper on separate stabilizing frame
- Two aluminum side rails / push handles
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

#### Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

#### Drawer Configuration:

- (2) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

#### Assembled Cart Dimensions

- 36.5" H x 32" W x 22" D  
(92.7 cm H x 81.3 W x 55.9 D)
- Cart net weight: 115 lbs
- Approximate shipping weight: 140 lbs



**Color Shown is Teal  
Additional Colors Available**