



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

**Aluminum Universal Line Medium 5 Drawer Procedure Cart  
 with Key Lock  
 #AL808K5**

Harloff's 20 years of innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.

**Specifications**

**Construction:**

- Aluminum construction

**Features:**

- Internal pullout shelf on right side
- Key Lock with two keys
- 26.0" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumpers
- Extended caster footprint for stability when drawers are extended
- 5" (12.7 cm) ball-bearing premium casters, all full swivel- one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered wheels-on and mostly assembled
- Pre-drilled threaded holes for standard accessories
- Latex free
- Durable, hard-baked powder coat finish available in 16 colors
- Twelve year warranty

**Drawer Configuration:**

- (3) – 3.25" (8.3 cm)
- (1) – 6.5" (16.5 cm)
- (1) – 9.75" (24.8 cm)

**Internal Drawer Dimensions:**

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

**Assembled Cart Dimensions:**

- 41" H x 31" W x 24.75" D (104.1 cm x 78.7 cm x 62.9 cm)
- Cart weight: 105 pounds (47.7 kg)
- Approximate shipping weight: 135 pounds (62 kg)



**Sixteen available color choices:**

**\*Antimicrobial Finish Available**

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black