

Universal Line Wide 8 Drawer Cart, Standard Package with Key Lock #AL8W10K8

Harloff's value treatment carts offer many of the same features our customers love in our large carts at an economical price point. This economy cart has become our fastest growing small treatment cart.



Construction:

- All aluminum, light-weight construction

Features:

- Extended width single column cart
- Light-weight aluminum construction
- Internal pullout shelf on right side
- Key locking cart with 2 keys
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumper
- Extended caster footprint for stability when drawers are extended
- 5" ball-bearing premium casters, all full swivel- one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled and ready to use
- Pre-drilled treaded holes for standard accessories
- Latex free
- Durable, hard-baked powder coat finish available in 17 colors
- Twelve year manufacturer's warranty

Internal Drawer Dimensions:

- 30" W x 17" D

Drawer Configuration:

- (7) – 3.25" multi-purpose storage drawer
- (1) – 9.75" multi-purpose storage drawer

Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Assembled Cart Dimensions:

- 48.5" H x 38.25" W x 25" D
- Writing surface: 44"
- Cart weight: 135 pounds
- Approx. shipping weight: 170 pounds