



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

## Universal Line Wide Medication Cart with Key Lock, Standard Package (MEDICINE ON TIME® Compatible)

### #AL8W10K4MOT

Harloff's Universal Line of medication carts were designed to accommodate all types of medication packaging. Its modern appearance incorporates industry-leading standard features and offers multiple locking options.

### Specifications

#### Construction:

- All Aluminum construction

#### Features:

- Extended width single column medication cart with up to 250 MOT® card capacity
- Locking narcotics storage included
- Internal pullout shelf on right side
- Key locking cart with 4 keys per lock
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumper
- Extended caster footprint for stability when drawers are extended
- 5" ball-bearing premium casters, all full swivel, 2 with toe locking brake & 1 directional
- Ball –bearing full extension drawer slides
- Cart is delivered fully assembled and ready to use
- Pre-drilled hole locations for standard accessories
- Latex free
- Durable, hard-baked power coat finish available in 17 colors
- Twelve year manufacturer's warranty

#### Drawer Configuration:

- 1 ea – 3.25" multipurpose storage drawer
- 1 ea – 10" drawer with one aluminum riveted divider and locking narcotics box.
- 2 ea – 10" drawer with two aluminum riveted dividers.

#### Internal Drawer Dimensions:

- 30" W x 17" D

#### Assembled Cart Dimensions:

- 48.5" H x 38.25" W x 25" D
- Writing Surface: 44"
- Cart weight: 125 pounds
- Approximate shipping weight: 160 pounds



#### Seventeen available color choices:

(may be viewed at [www.harloff.com](http://www.harloff.com))

Beige, Burgundy, Light Blue, Navy, Forest Green, Light Gray, Mauve, Red, Teal, Yellow, White, Hammer Black, Hammer Red, Hammertone Blue, Hammertone Gray, Hammertone Green, and Sand

## Universal Line Wide 4 Drawer Medication Cart with Key Lock for DISPILL® Packaging #AL8W10K4DSP

Harloff's Universal Line of medication carts were designed to accommodate all types of medication packaging. Its modern appearance incorporates industry-leading standard features and offers multiple locking options.



### Construction:

- Aluminum, light-weight cabinet construction with steel drawers

### Features:

- Extended width single column medication cart with 100 to 150 DISPILL® card capacity
- Locking narcotics storage included
- Internal pullout shelf on right side
- Key locking cart with 4 keys per lock
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumper
- Extended caster footprint for stability when drawers are extended
- 5" ball-bearing premium casters, all full swivel- one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled and ready to use
- Pre-drilled hole locations for standard accessories
- Latex free
- Durable, hard-baked powder coat finish available in 17 colors
- Twelve year manufacturer's warranty

### Internal Drawer Dimensions:

- 30" W x 17" D

### Drawer Configuration:

- (2) – 3.25" multipurpose storage drawer
- (1) – 13" drawer with locking narcotics box and rail system to hold DISPILL® packaging
- (1) – 13" drawer with rail system to hold DISPILL® packaging

### Assembled Cart Dimensions:

- 48.5" H x 38.25" W x 25" D
- Writing Surface: 44"
- Cart weight: 140 pounds
- Approximate shipping weight: 175 pounds

### Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

**Aluminum Universal Line Wide 4 Drawer Medication Cart for  
 DISPILL® Packaging with CompX® Electronic Lock  
 #AL8W10EKC4DSP**

Harloff’s Universal Line of medication carts were designed to accommodate all types of medication packaging. Its modern appearance incorporates industry-leading standard features and offers multiple locking options.

**Specifications**

**Construction:**

- Aluminum, light-weight cabinet construction with 13” steel drawers

**Features:**

- Extended width single column medication cart with 100 to 150 DISPILL® card capacity
- Locking narcotics storage included
- Internal pullout shelf on right side
- CompX® brand electronic keypad lock on cabinet with key lock override
- Access rights may be assigned for up to 250 users per cart
- Drawers automatically re-lock after access (pre-set programmed time)
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumper
- Extended caster footprint for stability when drawers are open
- 5” ball-bearing premium casters, all full swivel- one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled and ready to use
- Pre-drilled hole locations for standard accessories
- Latex free
- Durable powder coat finish available in 17 colors
- Twelve year manufacturer’s warranty



**Internal Drawer Dimensions:**

- 30” W x 17” D

**Assembled Cart Dimensions:**

- 48.5” H x 38.25” W x 25” D
- Writing Surface: 44”
- Cart weight: 145 pounds
- Approximate shipping weight: 180 pounds

**Seventeen available color choices:**

(may be viewed at [www.harloff.com](http://www.harloff.com))

Beige, Burgundy, Light Blue, Navy, Forest Green, Light Gray, Mauve, Red, Teal, Yellow, White, Hammer Black, Hammer Red, Hammertone Blue, Hammertone Gray, Hammertone Green, and Sand

**Additional Option: 41052 – Lockview® Software**

- USB dongle with RJ-11 cable.
- Allows user to interact with electronic lock through a laptop or PC.
- Software enables mass programming, view of audit trail (who accessed the cart and when), changes in auto re-lock timing, etc.

**Drawer Configuration:**

- (2) – 3.25” multipurpose storage drawer
- (1) – 13” drawer with locking narcotics box and rail system to hold DISPILL® packaging
- (1) – 13” drawer with rail system to hold DISPILL® packaging