

Pediatric / Model: M3030B9-PED-EMG

Modular 30 Line, Nine Drawer Pediatric Cart, EMG. PKG



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Sturdy uni-body steel Construction

Features:

- Tall cabinet with separate breakaway locks on each drawer
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cart is designed to be compatible with popular pediatric color coding systems
- Stabilizing pontoon bumpers
- 5" ball-bearing casters, all full swivel - one tracking, two with brake
- Antimicrobial plastic top, removable
- Flip down drawer handles
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Cart is delivered assembled
- Twelve year warranty on manufactured parts

Accessory Package Includes:

- IV Pole, Shipped Loose
- Aluminum O2 Holder, Shipped inside bottom drawer of cart
- Cardiac Board, Shipped Installed
- Side Rail Push Handle, Shipped Installed
- Utility Hooks Qty. 4, Shipped Installed

Drawer Configuration:

- 8 – 3" H (7.6 cm)
- 1 – 6" H (15.2 cm)

Internal Drawer Dimensions:

- 24.5" W x 17.25" D (62.3 cm W x 43.8 cm D)

Cart Dimensions:

- Overall Cart dims: 38.75" W x 22.25" D x 68.25" H (98.42 cm W x 56.51 cm D x 173.35 cm H)
- Empty Cart Weight: 155 lbs. (70.3 kg)
- Approximated Shipping Weight: 186 lbs. (84.3 kg)
- Approximated Shipping Dims: 37"L x 26"W x 56"H

Pediatric / Model: M3030B9-PED

Modular 30 Line, Nine Drawer Pediatric Cart, STD. PKG



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Sturdy uni-body steel Construction

Features:

- Tall cabinet with separate breakaway locks on each drawer
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cart is designed to be compatible with popular pediatric color coding systems
- Stabilizing pontoon bumpers
- 5" ball-bearing casters, all full swivel - one tracking, two with brake
- Antimicrobial plastic top, removable
- Painted steel top rail / push handle
- Flip down drawer handles
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Cart is delivered assembled
- Twelve year warranty on manufactured parts

Accessory Package Includes:

- IV pole (6804010)
- Aluminum O2 holder (AL680403)
- Cardiac board (680406)
- Utility hooks, (4) (680402)
- Customer installs accessories

Drawer Configuration:

- 8 - 3" H (7.6cm)
- 1 - 6" H (15.2cm)

Internal Drawer Dimensions:

- 24.5" W x 17.25" D (62.3cm W x 43.8 cm D)

Cart Dimensions:

- Overall Cart dims: 35.25" W x 22.25" D x 45.25" H (89.53 cm W x 56.51 cm D x 114.93 cm H)
- Empty Cart Weight: 151 lbs. (68.4 kg)
- Approximated Shipping Weight: 171 lbs. (77.5 kg)
- Approximated Shipping Dims: 37"L x 26"W x 56"H