



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Classic Line Medication Cart for Multi-Dose Packaging with Key Lock #MD30-4K

Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Key lock with two keys
- Vinyl-padded full wrap-around steel bumper
- Five-inch sealed ball-bearing premium casters, all full swivel – one with brake, one directional
- Internal locking narcotics box in 1/3 of bottom drawer
- Pull-out side shelf
- Aluminum push handle/side rail
- Removable plastic top
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled
- Certified "Latex Free"
- Five year warranty
- Compatible with most accessories
- Pre-threaded accessory mounting holes
- Hard-baked powder coat finish

- Available in Specialty Package. (MD30-4KSP)
 "SP" package includes sharps and waste container with brackets.

Drawer Configuration:

- 3 – 4" Drawers with plastic drawer dividers (10 spaces per drawer)
- 1 – 12" Drawer with 1 adjustable divider and locking narcotics box

Assembled Cart Dimensions:

- 38.5" H x 32" W x 22" D
 (97.8 cm H x 81.3 cm W x 55.9 cm D)
- Cart weight: 130 pounds (59 kg)
- Shipping weight: 155 pounds (61 kg)



Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
 Light Blue, Navy, Light Gray,
 Mauve, Red, Teal, Yellow,
 Hammer Tone Blue, Hammer Tone Gray
 and Hammer Tone Green