



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

3/4 Height Preconfigured Double Wide Catheter Procedure Cart with Tambour Door

#MS6140-CATH3

The MS6140-CATH3 is a ¾ height full shell, double wide Catheter Procedure Cart with tambour door, specifically designed for use in a cath lab to support the catheterization procedure. At five feet tall, the cart is 110 pounds lighter and over a foot shorter than the MS8140CATH

The MS6140-CATH3 comes preconfigured as follows: left column, top to bottom 1 catheter slide shelf with 3 slides; right column, top to bottom, 3ea 4" trays, 3ea 8" trays



Assembled Cart Dimensions:

- 60.5" H x 42" W x 27.9" D
(153.7 cm H x 106.7 cm W x 70.8 cm D)
- Cart net weight: 240 pounds (108 kg)
- Approximate shipping weight: 370 pounds (167 kg)

Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- The significant amount of catheter sizes and packaging techniques has made organizing and managing this critical area an extremely difficult task. Harloff offers a solution with a high density Catheter Storage System, which provides quick access and retrieval of catheters and related supplies.
- Equipped with adjustable Cath Slide Shelf, consisting of seven full extension catheter slides, each slide has five hooks, and each hook can hold up to 15 catheter packages.
- Each Cath Slide features a pair of ball bearing slides to withstand the pull force of removing catheters from hooks; full extension brings the catheters into full view for quick access and visual inventory of remaining stock.
- Cath Slide Shelf can be mounted in top position or 13 inches lower for storage above.
- Stores long stent boxes, balloons, guidewires, and other supplies efficiently for easy access with our specially-engineered catheter.
- The maximum package length is 55 inches.

Standard Features:

- Catheter Slide Shelf, featuring three cath slides and five hooks per slide, can be mounted in a lower position to create storage above for over-sized items. Cath Slide Shelf consists of a pair of ball bearing drawer slides for extra strength and durability
- Four inch modular trays: Total of three trays that are 23.50 x 15.75 x 4 inches
- Eight inch modular trays: Total of three trays that are 23.50 x 15.75 x 8 inches
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with cart side
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Double Column $\frac{3}{4}$ Height Mobile Medical Storage Cabinet with Key Locking Roll-up Tambour Door #MS6140CATH-4

The MS4040 $\frac{3}{4}$ height double column medical storage cabinet secured by a key locking roll-up tambour door. Utilizing the most innovative modular tray system on the market, this mobile cabinet can be configured for use in Arthroscopy, Cardiology, Catheterization, Cystoscopy, Endovascular Labs, ENT, Interventional Radiology, Laparoscopy, Ophthalmology, Orthopedic, Neurology, Urology and Vascular.

The MS6140-CATH4 comes preconfigured as follows: left column, top to bottom, 1 CATH-3 catheter slide shelf with 3 slides; right column, top to bottom, 2 ea 3" wired baskets 2 ea 5" wired baskets 2 ea 8" wired basket



Color Shown is SAND
Additional Colors Available

Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- Maximize space and minimize costs using 2", 4" and 8" modular trays; or 3", 5" and 8" wired baskets.
- Increased storage, reducing capital purchase cost.
- Cath Slide Shelf can be mounted in top position or 13 inches lower for storage above.
- Stores long stent boxes, balloons, guidewires, and other supplies efficiently for easy access
- The maximum package length is 55 inches.

Standard Features:

- Catheter Slide Shelf, featuring three cath slides and five hooks per slide
- Key locking roll-up tambour door
- Painted steel top with integrated push handles
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand
- Twelve year warranty against manufactures defect

Assembled Cart Dimensions:

- 60.5" H x 42" W x 27.9" D
(153.7 cm H x 106.7 cm W x 70.8 cm D)
- Cart net weight: 240 pounds (108 kg)
- Approximate shipping weight: 370 pounds (167 kg)