



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Standard Line 480 Punch Card Medication Cart, Standard Package
 #SL480PC**

Specifications



Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Dual column medication cart with 480 punch card capacity, single width narcotics drawer
- BEST® brand lock on cabinet and Illinois lock on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Offered in specialty package (#SL480PCSP) which includes a #68630-P1 drawer tray with adjustable dividers, #680551 eye drop tray, #AL680435 aluminum waste container and #684810 locking sharps container

Drawer Configuration:

Left Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 3 ea - 9.75" punch card drawers with adjustable dividers

Right Column Drawers:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

Internal Drawer Dimensions:

- Left Side: 17" D x 20.5" W (43.2 cm x 52.07 cm)
- Right Side: 17" D x 6.25" W (43.2 cm x 15.9 cm)

Thirteen available color choices:

Sand, Gloss Beige, Burgundy,
 Light Blue, Navy, Light Gray,
 Mauve, Red, Teal, Yellow,
 Hammer Tone Blue, Hammer Tone Gray
 and Hammer Tone Green

Assembled Cart Dimensions:

- 45.75" H x 36" W x 23.5" D
 (116.2 cm x 91.4 cm x 59.7 cm)
- Writing Surface: 44" H (111.8 cm)
- Cart weight: 170 lbs (77.3 kg)
- Approx shipping weight: 215 lbs (97.7 kg)