

Standard Line Punch Card Medication Cart, Standard Package #SL550PC



Construction:

- Painted steel cabinet and aluminum drawers

Features:

- Single column medication cart with up to 550 punch card capacity, single width internal narcotics box
- BEST® brand lock with 4 keys on cabinet.
- Illinois lock with 4 keys on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- 12 Year Manufacturer warranty

Internal Drawer Dimensions:

- 17" D x 29.5" W

Drawer Configuration:

- 1 ea - 4" multipurpose storage drawer
- 1 ea - 9.75" punch card drawer with dividers and locking narcotics box
- 2 ea - 9.75" punch card drawers with dividers

Assembled Cart Dimensions:

- 45.75" H x 36" W x 23.5" D
(116.2 cm x 91.4 cm x 59.7 cm)
- Writing Surface: 44" (111.8)
- Cart weight: 170 lbs (77.3 kg)
- Approx shipping weight: 210 lbs (95.5 kg)

Seventeen Available Colors:

Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand