



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Wooden Laminate 108 Unit-Dose Medication Cart, Standard Package #WL108BOX-SO

This wood laminate Unit-Dose cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.

Specifications



Construction:

- Aluminum base and internal frame
- MDF panels covered with a special heat-pressed thermal foil laminate with a rich wood appearance

Features:

- Single column medication cart with 108 unit-dose box capacity
- BEST® brand lock on cabinet
- Internal drawer dividers included in 6- 3.25" drawers
- Internal locking narcotics box in 1/2 of bottom 9.75" drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- 4" ball-bearing premium rubber casters all full swivel
- Bottom skirting conceals casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture

Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 6 ea – 3.25" drawers with internal dividers
- 1 ea – 9.75" punch card drawer with internal locking narcotics box in 1/2 drawer



Cherry Mahogany (CM)

WL108BOX-CM



Solara Oak (SO)

WL108BOX-SO

Color Options:

Color shown in Solara Oak with a gray granite finish on top. Also available in Cherry Mahogany #WL108Box-CM

Assembled Cart Dimensions:

- 44" H x 22.5" W x 22.5" D
- Writing Surface: 43.5" (110.5 cm)
- Cart weight: 175 lbs (79.5 kg)
- Approx shipping weight: 215 lbs (97.7 kg)

Additional laminate color options available with minimum cabinet purchase and extended lead time.