



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Wooden Laminate 240 Capacity Punch Card Medication Cart,
 Standard Package
 #WL240PC-CM**



Specifications

This wood laminate punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system that meets the demands of Long Term Care.

Construction:

- Aluminum base and internal frame
- MDF panels covered with a special heat-pressed thermal foil laminate with a rich wood appearance

Features:

- Single column medication cart with 240 punch card capacity
- BEST® brand lock on cabinet
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- 4" ball-bearing premium rubber casters all full swivel
- Bottom skirting conceals casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture

Narcotics Storage Option Available:

- Narcotics storage available.
WL240PCN-CM This provides an internal locking narcotics box with an Illinois lock in a 9.75" drawers.

Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 3 ea - 9.75" punch card drawer with adjustable dividers

Assembled Cart Dimensions:

- 44" H x 21.7" W x 22.5" D
- Writing Surface: 43" (109.2 cm)
- Cart weight: 160 lbs (72.7 kg)
- Approx shipping weight: 205 lbs (93.2 kg)

Available Color Choices and Options:



Cherry Mahogany (CM)

WL240PC-CM



Solara Oak (SO)

WL240PC-SO

- Cart above is shown in Cherry Mahogany finish
- The standard top color is Grey Granite
- Additional laminate color options available but require minimum purchase and extended lead times