



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

**Wooden Laminate 294 Capacity Unit Dose Medication Cart,  
 Standard Package  
 #WL294BOX-CM**

**This wood look mobile unit dose cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.**

**Specifications**



**Construction:**

- Aluminum base and internal frame
- MDF panels covered with a special heat-pressed thermal foil laminate with a rich wood appearance

**Features:**


- Single column medication cart with 294 unit dose capacity
- BEST® brand lock on cabinet and narcotics box located in bottom drawer
- Internal drawer dividers included in 6 – 3.25” drawers
- Adjustable row divider in bottom drawer for liquid storage
- Ball-bearing, full extension drawer slides
- Painted steel pull-out shelf
- 4” ball-bearing premium rubber casters all full swivel
- Bottom skirting conceals casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture.

**Narrow Cabinet available, #WL108BOX-CM or WL108BOX-SO**



Cherry Mahogany (CM)

WL480PC-CM



Solara Oak (SO)

WL480PC-SO

**Color Options:**

Color shown in Cherry Mahogany with a gray granite finish on top. Also available in Solara Oak #WL480PC-SO

Additional laminate color options available with minimum cabinet purchase and extended lead time.

**Drawer Configuration:**

- 1 ea – 4” multipurpose storage drawer
- 6 ea – 3.25” drawers with drawer dividers
- 1 ea – 9.75” drawer with locking narcotics box and adjustable divider

**Assembled Cart Dimensions:**

- 45” H x 35.5” W x 22.5” D
- Writing Surface: 44” (111.8cm)
- Cart weight: 230 lbs (107.8 kg)
- Approx shipping weight: 300 lbs (136 kg)