



650 Ford Street  
 Colorado Springs, CO 80915  
 Toll:(800)433-4064 \* Ph: (719)637-0300  
 Fax:(719)597-8273 \* www.harloff.com

**Wooden Laminate Treatment Cart, Standard Package  
 #WL6K-CM**



**Specifications**

**This wood laminate treatment cart provides an attractive alternative to traditional treatment carts. Designed to look like a piece of furniture with the internal design and locking system that meets the demands of Long Term Care.**

**Construction:**

- Aluminum base and internal frame
- MDF panels covered with a special heat-pressed thermal foil laminate with a rich wood appearance

**Features:**

- Single column treatment cart
- BEST® brand lock on cabinet
- Ball-bearing, full extension drawer slides
- 4" ball-bearing premium rubber casters all full swivel
- Bottom skirting conceals casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture

**Drawer Configuration:**

- 1 ea – 4" storage drawer
- 3 ea- 3.25" storage drawers
- 2 ea - 9.75" storage drawer

**Assembled Cart Dimensions:**

- 44" H x 21.7" W x 22.5" D (111.8cm x 57.8cm x 57.2cm)
- Writing Surface: 43.5" (110.5 cm)
- Cart weight: 175 lbs (79.5 kg)
- Approx shipping weight: 215 lbs (97.7 kg)

**Available Color Choices and Options:**



Cherry Mahogany (CM)



Solara Oak (SO)

WL6K-CM

WL6K-SO

- Cart above is shown in Cherry Mahogany finish
- The standard top color is Grey Granite
- Additional laminate color options available but require minimum purchase and extended lead times