



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Wood Vinyl 360 Capacity Punch Card Medication Cart, Standard Package #WV360PC-CM

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.



Specifications

Construction:

- Vinyl coated aluminum
- Top is MDF covered with a heat-pressed thermal foil laminate

Features:

- Light weight vinyl coated aluminum
- Single column medication cart with 360 punch card capacity
- Key locking cabinet with locking narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- Pull-out shelf
- 5" ball-bearing premium rubber casters, full swivel, 1 with toe locking brake and 1 directional
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture

Drawer Configuration:

- 1 ea – 3" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers



Cherry Mahogany (CM)

WV360PC-CM



Buffalo Oak (OAK)

WV360PC-OAK

Color Options:

The WV360PC-CM is shown above in Cherry Mahogany finish with a grey granite finish on the top. Additional vinyl color options available with minimum cabinet purchase and extended lead time.

Assembled Cart Dimensions:

- 44.5" H x 30" W x 23.5" D
- Writing Surface: 44"
- Cart weight: 112 lbs
- Approx shipping weight: 142 lbs