



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Wood Vinyl 540 Capacity Punch Card Medication Cart, Standard Package #WV540PC-CM

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.

Specifications



Construction:

- Vinyl coated aluminum

Features:

- Light weight vinyl coated aluminum
- Single column medication cart with 540 punch card capacity
- Key locking cabinet and narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- Pull-out shelf
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture

**Additional Wood Look Option,
 #WL480PC-CM Wood Laminate with Bottom Skirting**



Cherry Mahogany (CM)

WV540PC-CM



Buffalo Oak (OAK)

WV540PC-OAK

Color Options:

The WV540PC-CM is shown above in Cherry Mahogany finish with a grey granite finish on the top. Additional vinyl color options available with minimum cabinet purchase and extended lead time.

Drawer Configuration:

- 1 ea – 3" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers

Assembled Cart Dimensions:

- 45" H x 34.5" W x 24.5" D
- Writing Surface: 44" (111.8)
- Cart weight: 127 lbs (57.7 kg)
- Approx shipping weight: 157 lbs (71.2 kg)



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

Wood Vinyl 540 Capacity Electronic Locking Punch Card Medication Cart, Standard Package #WV540PCEKC-CM

This vinyl coated aluminum punch card cart provides an attractive alternative to traditional medication carts. Designed to look like a piece of furniture with the internal design and locking system to meet the demands of Long Term Care use.



Specifications

Construction:

- Vinyl coated aluminum

Features:

- Light weight vinyl coated aluminum medication cart with 540 punch card capacity
- Electronic locking cart with keypad access and key lock override
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart
- May add, edit, or delete user and user groups
- Battery operated with six D-cell Batteries; life 1 year minimum at 25-50 uses/day
- Lock has an LED low battery indicator
- Narcotics box located in top drawer
- Adjustable row dividers for punch card and liquid storage
- Ball-bearing, full extension drawer slides
- Pull-out shelf
- 5" ball-bearing premium rubber casters
- Cart is washable with soap and water
- This cart is designed to look like a piece of furniture



Cherry Mahogany (CM)

WV540PCEKC-CM



Buffalo Oak (OAK)

WV540PCEKC-OAK

Color Options:

The WV540PCEKC-CM is shown above in Cherry Mahogany finish with a grey granite finish on the top. Additional vinyl color options available with minimum cabinet purchase and extended lead time.

Assembled Cart Dimensions:

- 45" H x 34.5" W x 24.5" D
- Writing Surface: 44" (111.8)
- Cart weight: 127 lbs (57.7 kg)
- Approx shipping weight: 157 lbs (71.2 kg)

Additional Wood Look Option, #WL480PC-CM Wood Laminate with Bottom Skirting

Lockview Software Available:
 #41052 Includes audit software and dongle



Drawer Configuration:

- 1 ea – 3" multipurpose storage drawer
- 1 ea – 9.75" punch card drawer with locking narcotics and adjustable dividers
- 2 ea - 9.75" punch card drawer with adjustable dividers