

A. Model 1362



Common Features

POWER PARALLEL BARS —
MODELS 1361, 1362, 1363, 1364:



- Motorized height adjustment from 27" to 37".
- Convenient hand control adjusts bars.
- 1 1/8" diameter stainless steel handrails.
- 34" wide wood platform is covered in textured gray styrene.
- Tapered ends for wheelchair access.
- Removable Abduction Board.
- Manual crank override.
- Steel uprights are powder coated in neutral cream finish.
- 30 1/2" clearance between uprights.
- Hospital grade plug.
- 110 V, 60Hz, 3 Amps. U.L Listed.

A. Models 1362, 1363 Power Height and Width Parallel Bars

SPECIAL FEATURES:

- Motorized width adjustment of handrails from 13" to 25".
- Model 1363 15' long has (6) uprights.

MODEL	LENGTH	WIDTH	HEIGHT
1362	10'	13"- 25"	27"- 37"
1363 (not shown)	15'	13"- 25"	27"- 37"

Models 1361, 1364 Power Height Parallel Bars

SPECIAL FEATURES:

- Can be used by adults and children.
- Width adjusts manually from 15" to 27".
- Includes brackets to lower handrails to child height.*
- Model 1364 15' long has (6) uprights.

MODEL	LENGTH	WIDTH	HEIGHT
1361 (not shown)	10'	15"- 27"	27"- 37"*
1364 (not shown)	15'	15"- 27"	27"- 37"*

*child height is 17"- 27"

B. Models 1365, 1366 Crank Height Platform Bars

FEATURES:



- Hand crank height adjustment of each stainless steel handrail 26" to 39".
- 34" wide wood platform is covered in textured gray styrene.
- Tapered ends for wheelchair access.
- Hand controls allow easy adjustment of width from 16 1/2" to 24 1/2".
- Ruled height and width adjustments.
- Model 1366 15' long has (6) uprights.
- Includes removable Abduction Board.

MODEL	LENGTH	WIDTH	HEIGHT
1365	10'	16 1/2" - 24 1/2"	26" - 39"
1366*(not shown)	15'	16 1/2" - 24 1/2"	26" - 39"

*1366 not Quick Ship

B. Model 1365

