



A. Model 1390

Common Features

PLATFORM MOUNTED PARALLEL BARS—
MODELS 1300, 1384, 1390, 1391, 1392:



- Satin-finish hardwood platform with tapered hardwood ends for easy wheelchair access.
- (2) “anti-slip” treads on each end.
- 1½” diameter one-piece stainless steel handrails.
- Heavy gauge black powder coated steel uprights and fittings.
- Each upright telescopes up in 1½” increments and locks into (10) height positions with fail-safe ball-tip locking pin.
- 3mm Top edge for greater durability.

OPTION

FOR PLATFORM PARALLEL BARS

- *62 – Abduction Board
- *63 – Non-slip Matting

*FACTORY INSTALLED

A. Models 1390, 1391, 1392 Height & Width Adjustable Parallel Bars

SPECIAL FEATURES:

- Width adjusts from 15” to 28” with ergonomic control knobs on each upright.
- 30” clearance between uprights.

MODEL	LENGTH	WIDTH*	HEIGHT
1390	10'	15"- 28"	29"- 42"
1391 (not shown)	7'	15"- 28"	29"- 42"
1392 (not shown)	12'	15"- 28"	29"- 42"

*platform is 38”wide

Model 1300 (not shown) Height Adjustable Parallel Bars

MODEL	LENGTH	WIDTH*	HEIGHT
1300	10'	26"	28"- 41"

*platform is 38”wide

B. Model 1384 Combination Adult-Child Parallel Bars

FEATURES:

- Unit includes adjustable height Adult Bar and adjustable height Child’s Bar.
- Adult bar adjusts from 29” to 42”H and Child’s bar adjusts from 20” to 32” high.
- Each upright telescopes up/down in 1½” increments and locks into position with fail-safe ball-tip locking pin.
- Both Bars adjust in width from 15” to 28” with ergonomic control knobs.
- 10 ft long Satin-finish hardwood platform with tapered ends for easy wheelchair access.
- 30” fixed clearance between black powder-coated steel uprights.

MODEL	LENGTH	WIDTH*	HEIGHT
1384	10'	15"-28"	29"-42"



B. Model 1384