



A. MODEL 1427-57-121

## COMMON FEATURES

MAT PLATFORMS —

SERIES 1427, 1428, 1429:



- Extremely stable lift system elevates from 20" wheelchair height to 30" high.
- "Hands-free" foot control.
- Seamless nylon-reinforced vinyl is wear and tear resistant.
- 2" High-density urethane foam vinyl top.
- Steel frame and 4 corner legs have durable black powder coat finish.
- Legs have levelers for uneven floors.
- 4-Leg "open" design is patient sling-lift accessible on all four sides.
- Hospital grade plug.

**COLORS:** (SEE PAGE 24)

**Upholstery** – 113 Blue, 119 Gray, 121 Green, 122 Black

You must specify upholstery color with order.

### A. SERIES 1427 BARIATRIC MAT PLATFORMS



**SPECIAL FEATURES:**

- Electrically activated hydraulic cylinders are built into each (4) corner steel leg.
- Motorized lift system can securely lift 750 lbs. from 20" wheelchair height to 30" high.
- "Hands-free" foot control operation.

MODEL	WIDTH	LENGTH	HEIGHT
1427-46*	4'	6'	20"-30"
1427-47*	4'	7'	20"-30"
1427-57*	5'	7'	20"-30"
1427-68*	6'	8'	20"-30"

### B. SERIES 1428 ELECTRIC HI-LO MAT PLATFORM WITH POWERBACK



**SPECIAL FEATURES:**

- Electrically activated hydraulic cylinders built into each of 4 corner legs with 750 lbs. patient weight capacity.
- Wheelchair Accessible - Elevates smoothly from 21" wheelchair height to 31" H.
- Electric backrest can lift and elevate the back of an immobile patient. Backrest lift capacity is 400 lbs. Backrest adjusts from 1° to 70°.
- "Hands free" foot control operation for backrest and height adjustment.
- Ready for Use - Ships fully assembled ready for use. Just Plug it in and Good to Go!

MODEL	WIDTH	LENGTH	HEIGHT
1428-47	4'	7'	21"-31"
1428-57	5'	7'	21"-31"

B. MODEL 1428-47-119



Series 1428  
Electric  
Backrest



C. MODEL 1429-68-113

### C. SERIES 1429 BARIATRIC 1100 LB. CAPACITY HEAVY-DUTY MAT PLATFORM



**SPECIAL FEATURES:**

- Electrically activated hydraulic cylinders are built into each (4) corner steel leg.
- Motorized lift system can securely lift 1100 lbs. from 20" wheelchair height to 30" high.
- Static load tested to 4400 lbs.
- Seamless nylon reinforced vinyl is tear resistant.
- Plywood reinforced subtop for added strength.

MODEL	WIDTH	LENGTH	HEIGHT
1429-57	5'	7'	20"-30"
1429-68	6'	8'	20"-30"