

A. Model 1570
Open setup



Closed setup



A. Model 1570
Patented 3-in-1
Training Staircase



FEATURES:

- U.S. Patent No. D.428,629.
- Unique training staircase device allows the therapist to train patients on (3) progressive step levels (4" – 6" – 8").
- (6) 4" high steps on one end.
- (4) 6" high steps on other end.
- (3) 8" high steps at front side.
- Steps are 36" wide; platform is 36"x 36".
- Handrail for 8" stair section is easily removable and locks into place.
- Features BriteStep™ 3mm edging
- Requires minimum of 9 ft. high ceiling.

MODEL	LENGTH	WIDTH	HEIGHT
1570	109"	56"*	24"

*includes 8" high steps

B. Model 1564
Bariatric ADA
Two-Sided Staircase



FEATURES:

- Natural Oak high pressure laminated reinforced plywood construction for 750 lb. patient weight capacity.
- (2) 6" high steps on one side and (3) 4" high steps on opposite side with BriteStep™ 3mm edging.
- Steps are 38" wide x 11" A.D.A. Depth
- Platform is 38" wide x 27" deep x 12" high.
- Solid laminated plywood side panels with ergonomic easy-to-grip hardwood handrails.
- Handrails have 12" ADA extension beyond bottom steps.
- Shipped ready-to-assemble.

MODEL	LENGTH	WIDTH	HEIGHT
1564	84"*	40"	12"

*Including 12" handrail extensions

C. Model S-601
Econo Two-Sided
Staircase



FEATURES:

- Two 6" steps and three 4" steps.
- Steps and platform have "anti-slip" safety treads.
- Platform is 30" W x 24" D x 12" H
- Compact size can fit in 8' ceilings.
- Folkstone Gray laminate finish.
- Ergonomic hardwood handrails.

MODEL	LENGTH	WIDTH	HEIGHT
S-601	52"	30"	12"

C. Model S-601



B. Model 1564

