

A. Model 2153 Economy Air-Lift Stool



- Single lever height Control Handle.

MODEL	DIAMETER	HEIGHT
2153	16"	20" - 27"
2153-S <i>New Product!</i>	16"	17" - 22"

B. Model 2157 Economy Air-Lift Stool



- 360° circular ring Control Handle.

MODEL	DIAMETER	HEIGHT
2157	16"	20" - 27"
2157-S <i>New Product!</i>	16"	17" - 22"

Common Features

"CHOICE LINE™" AIR-LIFT STOOLS
MODELS 2153, 2157:



- Hand operated air-spring pneumatic height adjustment.
- 16" diameter seat, 3½" thick foam.
- (5) 2" Twin wheel black nylon casters.
- Durable 5-leg, black composite base.

COLORS: (FOLD OUT FRONT COVER)

Upholstery – 17 standard vinyl colors.



OPTIONS FOR MODELS 2153, 2157

- *L – Debossed Logo
Only available in PROFORM colors.
A900 6" plate for top debossing.
Scan QR code for more info on Logo Stools.



6" Logo Stool Top

*FACTORY INSTALLED

SEE PAGE 83

C. Model 2151 Heavy-Duty Air-Lift Stool



- Heavy-duty Pneumatic Lift with "D" chrome Control Handle.
- 16" Diameter Seat, 3½" thick foam seat.
- 4" stroke cylinder.
- Durable Aluminum 5-Leg Base.
- (5) 2½" Twin wheel gray nylon casters.
- Heavy-duty 400 lb. weight capacity.
- Available in 17 standard vinyl colors.

MODEL	DIAMETER	HEIGHT
2151	16"	20"-24"

D. Model 2154 "Dual" Foot & Hand Operated Stool



- Pneumatic Stool has BOTH Foot & Hand Controls for foot or hand operation.
- 5" stroke cylinder.
- Includes 360° ring Hand Control & circular Metal Foot Control.
- Durable 5-leg Polished Aluminum Base .
- (5) 2" Twin wheel black nylon casters.
- Available in 17 standard vinyl colors.

MODEL	DIAMETER	HEIGHT
2154	16"	19"-24"

E. Model 2156-707 Economy Air-Lift Stool with Backrest



- Padded backrest, 10" high x 15" wide.
- Air-spring pneumatic height adjustment.
- 360° circular ring Control Handle.
- Durable 5-leg, black composite base.
- 16" diameter seat, 3½" thick foam.
- (5) 2" Twin wheel black nylon casters.
- Black vinyl upholstery.

MODEL	DIAMETER	HEIGHT
2156-707	16"	20" - 27"

C. Model 2151



D. Model 2154



E. Model 2156-707

