

A. Model 2410



A. Model 2410
Nested Footstools



FEATURES:

- Consists of (4) steps which “nest” together for easy space-saving storage.
- Nested footstools are useful for stepping, crawling or climbing activities.
- Durable oak laminate wood construction with non-slip matting.

MODEL	LENGTH	WIDTH	HEIGHT
2410 (closed)	20"	16"	8"
Stool Sizes:	20"	16"	8"
	18"	14"	6"
	14"	12"	4"
	12"	9"	2"

B. Model 1200



B. Model 1200
Ramp & Curb Set




FEATURES:

- Facilitates ambulation and patient mobility training on ramp and curbs.
- Durable Oak laminate wood construction.
- Anti-slip safety matting on top surface of ramp and curb set.
- Ramp elevates from ¼" to 8" high over its 44" incline.
- Flat platform is 44" x 44".
- Unit includes (3) curbs of 2", 4" and 6" height, which nest under ramp.
- Curbs can be set up to extend in any of 3 directions or curbs can be removed from unit and used separately.
- Features BriteStep™ 3mm Yellow edging on all curbs.

MODEL	LENGTH	WIDTH	HEIGHT
1200 (closed)	88"	44"	8"
Curb sizes	40"	40"	6"
	36"	36"	4"
	32"	32"	2"

Curbs can be set up to extend to right as shown above or to the left (below) or straight.



C. Model 1310



New Design!

C. Model 1310
Balance Activity Platform



FEATURES:

- Adjustable height 3' stainless steel handrails.
- Black high-pressure laminate with safety treads.

MODEL	WIDTH	DEPTH	HEIGHT
1310	36"	36"	28" - 41"

