




A. POSITIONING PILLOWS



B. TABLE PILLOWS



A. Positioning Pillows



FEATURES:

- Firm High-Density urethane foam with vinyl sewn covering.

COLORS:

Color for the Positioning Pillows is Medium Blue.

MODEL	TYPE	LENGTH	WIDTH	HEIGHT
30	Round Bolster	24"	6" Dia.	–
31	Round Bolster	24"	8" Dia.	–
33	Cube	12"	8"	3"
34	Cube	12"	8"	6"
36	20° Wedge	18"	24"	6"
37	20° Wedge	32"	20"	12"
38	Cube w/40° edge	32"	10"	12"
39	Cube w/40° edge	32"	20"	12"

B. Table Pillows




- Firm urethane foam with Resist-All™ vinyl upholstery.

COLORS:

Upholstery – 17 standard vinyl colors.
5 PVC-Free urethane colors.

MODEL	TYPE	LENGTH	WIDTH	HEIGHT
32	Full Size Pillow	22"	14"	3"
35	Small Pillow	12"	14"	3"

C. "ANTI-SLIP" BOLSTERS



C. Models 41, 42, 43 "Anti-Slip" Positioning Bolsters



- Firm High-Density urethane foam with Medium Blue vinyl cover.
- Bolsters have "anti-slip" bottom surface that stays in place on mats or tables.

MODEL	TYPE	LENGTH	WIDTH	HEIGHT
41	45° Wedge	16"	8"	8"
42	Half-round	17"	8"	7"
43	20° Wedge	32"	20"	12"

D. Model 3004 Positioning Wedge System



- Versatile wedge system for positioning of children or adolescents.
- Consists of (2) large and (2) small vinyl covered foam wedges with VELCRO® brand fasteners.
- Wedges assemble in a variety of configurations.
- Vinyl color is Medium Blue.

MODEL	LENGTH	WIDTH	HEIGHT
3004	31"	20"	17"
Consists of: (2) Small Wedges			
	14"	20"	8"
(2) Large Wedges			
	27"	20"	14"

D. Model 3004



E. Model 3003 Multi-Function Twin Wedge Set



- Two "twin" High-Density urethane foam Wedges covered in Medium Blue vinyl.
- Each 20" x 20" x 10"H Wedge has VELCRO® brand fasteners so Wedges can be attached back-to-back or put together as a Cube.
- Wedges can also be used separately for head, leg or knee support.

MODEL	WIDTH	LENGTH	HEIGHT
3003	20"	20"	10"

E. Model 3003





A. Model 1670



B. Model 1671



C. Model 1672

**A. Model 1670
Portable Triple Mirror**



- (3) high quality, distortion-free, full view glass mirrors with ANSI safety backing.
- Oak laminate frame and finished back.
- Mounted on 2" diameter swiveling casters.
- Viewing angles are easily changed and locked in desired position. Each mirrored panel is 24" wide.

MODEL	WIDTH	DEPTH	HEIGHT
1670	60"(open)	24"(open)	69"

**B. Model 1671
Portable Single Mirror**



- High quality, distortion-free glass mirror with ANSI "shatter-stop" safety backing.
- Oak laminate frame and base.
- Finished oak laminate back.
- Mounted on easy rolling 2" diameter swiveling casters.
- Mirror height provides full body image.

MODEL	WIDTH	DEPTH	HEIGHT
1671	24"	18"	70"

**C. Model 1672
Wall Mounted Mirror**



- High quality, full view, distortion-free glass mirror w/ANSI "shatter-stop" safety backing.
- Oak laminate frame.
- Equipped with mounting cleat for attachment to wall.
- Mirror frame measures 66" x 24". Can be mounted 6" above floor so that overall height is 72".

MODEL	WIDTH	DEPTH	HEIGHT
1672	24"	4"	66"



D. Model S-410
Front View



E. Model S-420
Reverse view

**D. Model S-410
Econo Mirror**



- Full-view glass mirror with ANSI safety backing.
- Gray laminate with (4) 2½" casters.

MODEL	WIDTH	DEPTH	HEIGHT
S-410	24"	18"	70"

**E. Model S-420
Econo Mirror Rack**



- Same as S-410 but with Pegboard back and (11) hooks.

MODEL	WIDTH	DEPTH	HEIGHT
S-420	24"	18"	70"

**Model S-420-100
Accessorized Mirror Rack**



MODEL	CONSISTS OF
S-420-100	(1) S-420 mirror, (1) 5581 Set of (16) Huger Weights

All weights shown must be ordered as accessories