

# SPECIALTY TREATMENT TABLES



**A. Model 4079**

## A. Model 4079 Modality Treatment Table



### FEATURES:

- Designed to save time and effort with modality equipment at your fingertips.
- Spacious upholstered top with 2" High-Density urethane foam.
- Oak laminate with (2) fully enclosed Door Cabinets with adjustable shelf.
- Mid-section of table has full extension 10"H x 23¾"W equipment compartment.
- Black laminate base has access holes on both ends and back for cord and plug.
- Choice of (9) Pro-Form™ vinyl colors.

MODEL	LENGTH	DEPTH	HEIGHT
4079	78"	30"	31"

## B. Model 7604-752 Portable Massage Table



### FEATURES:

- 2½" thick multi-layered semi-firm foam.
- Comfortable PVC-Free urethane upholstery.
- Precision welded yet lightweight aluminum frame.
- Push-button height adjustment.
- 600 lbs. weight capacity.
- Carrying bag with shoulder sling.
- Adjustable face cradle with ultra soft face cushion & arm sling.
- Spacious 73" long x 28" wide top.
- Height adjustable from 24" to 32".
- Color: 752 Indigo Blue Urethane.
- Weighs only 33 lbs.

MODEL	LENGTH	WIDTH	HEIGHT
7604-752	73"	28"	24" - 32"



**B. Model 7604-752**



**C. Model S-110**

## C. Model S-109, S-110 Econo-Line Table



### FEATURES:

- Oak laminate on all surfaces.
- Solid end panels with full-width storage shelf.
- Upholstered flat top has 2" urethane foam and choice of 706 Blue or 709 Gray vinyl.
- Ships ready-to-assemble.

MODEL	LENGTH	WIDTH	HEIGHT	UPH.COLOR
S-109	72"	27"	31"	706 Blue
S-110	72"	27"	31"	709 Gray

## OPTIONS

FOR S-109, S-110

- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small

SEE PAGE 4