

## Common Features

EXAM TABLES — MODELS 4416, 4412:



- Welded heavy-duty steel chassis with durable epoxy powder coated Gray Beige finish.
- One-piece vinyl upholstered top with High-Density urethane foam and contoured padded head section for patient comfort.
- Extra-wide 30" Top accommodates larger patients.
- Table height is comfortable 33" with toe recess on both sides of table.
- Spacious Footstep is 19½" W x 17" D x 9" H with anti-slip treaded surface and front supports for stability.
- Table is equipped with (2) front end Drawers and (3) right side Drawers. Molded polystyrene interiors for easy cleaning.
- Laminate Drawer fronts are ¾" thick with white 3mm edges for protection and durability.
- Heavy duty concealed Stirrups adjust in length and are laterally adjustable into 4 positions. Heel rests are powder coated for patient comfort.
- Slide-out 20" wide leg rest with Pad extends full 18".
- Pelvic Tilt with easy-to-use dual side controls.
- Plastic Irrigation Pan extends and retracts from legrest.
- Includes Paper Dispenser and Cutter.

**COLORS:** (FOLD OUT FRONT COVER)

**Upholstery** – 17 standard vinyl colors available.

### A. Model 4416 X-L Power-Back Exam Table



**Super Strong Steel Table for Bariatric and General Practice Use**

**SPECIAL FEATURES:**

- Backrest is motorized 0° to 75° to facilitate Bariatric or Elderly patient positioning.
- 600 lb. patient weight capacity.
- Electrical outlet in rear of Table with Hospital grade plug.
- 110V, 60Hz, 5 Amps.

MODEL	LENGTH	WIDTH	HEIGHT
4416	54" - 72"*	30"	33"

\*Extended

### B. Model 4412 Gas Spring Back Exam Table



**SPECIAL FEATURES:**

- Backrest is gas spring activated with control levers on both sides, adjusts from 0° to 75° and locks into any position.
- 450 lb. patient weight capacity.

MODEL	LENGTH	WIDTH	HEIGHT
4412	54" - 72"*	30"	33"

\*Extended

## OPTIONS

FOR MODELS 4412, 4416:

\*99L – Left Hand Table Model

\*FACTORY INSTALLED

SEE PAGE 4

A. Model 4416



Power Backrest is standard on 4416



Oversize 600 lb. rated steel footstep

B. Model 4412

