

SPECIALTY TREATMENT TABLES

Common Features

PROFESSIONAL TREATMENT TABLES —
MODELS 4524, 4541, 4543:



- 20" wide laminate drawer on nylon rollers with positive stops.
- 2" High-Density urethane foam box-style upholstered top.
- Resist-All™ vinyl upholstery finish.

COLORS:

Laminate – Natural Oak, Wild Cherry,
Folkstone Gray

Upholstery – 17 standard vinyl colors


(FOLD OUT FRONT COVER)



A. Model 4524

OPTIONS

FOR 4500 SERIES TABLES

- *00 – Special Table Height (N/A Model 4543)
- *04 – Table Safety Strap
- *06 – Backrest, adjustable 
- *07 – Nose Cutout
- *19 – Slide Shelf
- 23 – Paper Dispenser
- 25 – Paper Cutter
- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small
- *99L – Left hand table model

*FACTORY INSTALLED

SEE PAGE 4

A. Model 4524

SPECIAL FEATURES:

- H-Brace with laminate storage shelf.
- High pressure laminate legs, apron and stretchers.

MODEL	LENGTH	WIDTH	HEIGHT
4524	78"	30"	31"

B. Model 4541

SPECIAL FEATURES:

- H-Brace with laminate shelf and open shelf storage cabinet.
- High pressure laminate legs, apron and stretchers.

MODEL	LENGTH	WIDTH	HEIGHT
4541	78"	30"	31"

C. Model 4543

SPECIAL FEATURES:

- Fully enclosed all laminate storage cabinet with sliding laminate doors.
- Laminate interiors with divider and (1) adjustable shelf.
- Black laminated plywood base.

MODEL	LENGTH	WIDTH	HEIGHT
4543	78"	30"	31"



B. Model 4541

C. Model 4543

