


## Common Features

**HI-LO TABLES – MODELS 4701, 4702, 4703, 4704:**

- Hydraulic cylinders are built into each (4) Corner Telescoping Legs. 
- Black powder coated Legs with Levelers for uneven floors.
- Vinyl upholstered tops with 2" High-Density urethane foam.
- Includes Paper Dispenser and Cutter.
- 16"L legrest extends overall top length to 68".

**COLORS: (FOLD OUT FRONT COVER)**

**Upholstery – 9 Pro-Form™ vinyl colors.**



**A. Model 4703-722**

### A. Model 4703 Crank Hydraulic Changing/ Treatment Table



**SPECIAL FEATURES:**

- Crank hydraulic adjustment from 20" wheelchair height to 30" high.

MODEL	LENGTH	WIDTH	HEIGHT
4703	52" - 68"	30"	20" - 30"
4703-722 	52" - 68"	30"	20" - 30"

### B. Model 4702 Electric Hi-Lo Changing/ Treatment Table



**SPECIAL FEATURES:**

- Electrically activated hydraulic cylinders are built into each (4) Corner Telescoping Legs.
- Motorized table lowers to 20" wheelchair height. Elevates up to 30" high.
- "Hands-Free" foot control operation.

MODEL	LENGTH	WIDTH	HEIGHT
4702	52"-68"	30"	20" - 30"



**B. Model 4702**

### C. Model 4701 Crank Hydraulic Changing/ Treatment Table with Drawers



**SPECIAL FEATURES:**

- Crank hydraulic adjustment from 27"H to 37"H.
- Hydraulic cylinders are built into each (4) corner telescoping legs.

MODEL	LENGTH	WIDTH	HEIGHT
4701	52"-68"	30"	27" - 37"

### Model 4704 (not shown) Electric Hi-Lo Changing/ Treatment Table w/ Drawers



**SPECIAL FEATURES:**

- Electrically activated hydraulic cylinders are built into each (4) Corner Telescoping Legs.
- Motorized table lowers to 27". Elevates up to 37"H.
- "Hands-Free" foot control operation.

MODEL	LENGTH	WIDTH	HEIGHT
4704 (not shown)	52"-68"	30"	27" - 37"



**C. Model 4701**