

Common Features

POWER TABLES — MODELS 4717, 4719, 4723:



- Heavy-duty, motorized hydraulic Hi-Lo system with 10" lift stroke.
- Electrically activated hydraulic cylinders are built into (4) corner telescoping legs.
- Legs are black powder-coated and equipped with levelers for uneven floors.
- "Hands-Free" foot control operation.
- Vinyl upholstered tops with 2" High-Density urethane foam.
- 110V, 60Hz, 2.4 amps.

COLORS: (FOLD OUT FRONT COVER)

Upholstery – 9 Pro-Form™ vinyl colors.

A. Model 4723

ADA Hi-Lo Power Plinth Table



SPECIAL FEATURES:

- Bariatric 600 lb. lift capacity.
- Motorized table lowers to 20" for wheelchair transfer. Elevates up to 30" high.
- Sturdy black powder-coated steel frame with black laminate storage shelf.

MODEL	LENGTH	WIDTH	HEIGHT
4723	72"	30"	20"-30"

B. Model 4717

Hi-Lo Power Cabinet Treatment Table



SPECIAL FEATURES:

- Bariatric 500 lb. lift capacity.
- Motorized table elevates from 27" to 37" high.
- Fully enclosed oak laminate cabinet, finished inside and out, with sliding oak laminate doors.

MODEL	LENGTH	WIDTH	HEIGHT
4717	78"	30"	27"-37"

C. Model 4719

Hi-Lo Power Plinth Table with Shelf



SPECIAL FEATURES:

- Bariatric 600 lb. lift capacity.
- Motorized table elevates from 27" to 37" high.
- Oak laminate apron with black laminate storage shelf.

MODEL	LENGTH	WIDTH	HEIGHT
4719	78"	30"	27"-37"

OPTIONS

FOR MODELS 4717, 4719, 4723

- *15 – (1) Safety Rail (4717 only)
- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small

SEE PAGE 4



A. Model 4723
*Shown with option:
27 – Paper Dispenser & Cutter*



B. Model 4717



C. Model 4719
*Shown with option:
27 – Paper Dispenser & Cutter*