



**A. Model 4790 with Drop-Away "Open"**  
05 – Headrest is now standard.

### Common Features

ECHO-SCAN™ ECHOCARDIOGRAPH TABLES MODELS 4790, 4892:



- 3" High-Density urethane foam top.
- 10" wide x 12" deep drop-away section on durable piano hinge.
- Sliding latch activates drop-away section.

**COLORS:**

Upholstery – 17 standard vinyl colors  
*(FOLD OUT FRONT COVER)*

### A. Model 4790 Powermatic® Echo-Scan™ Table



**SPECIAL FEATURES:**

- Two direct-drive synchronized satin aluminum lift mechanisms.
- Power height adjustment from 20" wheelchair height to 29" high.
- "Hands-Free" foot control.
- Hospital grade plug.
- 110V, 60 Hz, 2.8 Amps. U.L. Listed.
- Includes 05 Headrest.

MODEL	LENGTH	WIDTH	HEIGHT
4790	76"	27"	20"-29"



**Model 4790 with Drop-Away "Closed"**  
*Shown with option: 074 - 4" Locking Casters (not intended for patient transport), 27 – Paper Dispenser & Cutter*



**B. Model 4892 with Drop-Away "Open"**

### B. Model 4892 Echo-Scan™ Table

**SPECIAL FEATURES:**

- Two laminated solid panel legs.
- Convenient laminate storage shelf.
- Folkstone Gray laminate finish.

MODEL	LENGTH	WIDTH	HEIGHT
4892	76"	27"	31"

### OPTIONS

FOR MODELS 4790, 4892

- \*04 – Table Safety Strap
- \*05 – Headrest, manually adjustable (4892 only)
- 27 – Paper Dispenser & Cutter Combo
- 36 – Wedge Pillow (see page 45)
- \*074 – (4) 4" locking casters (4790 only) (adds 4" to unit height)

\*FACTORY INSTALLED

SEE PAGE 4