

WEIGHT RACKS



A. Model 5570



B. Model 5565 Front view



Model 5565 Reverse view

All weights and dumbbells shown must be ordered as accessories

A. Model 5570 Weight Rack



- Almond coated pegboard on both front and back side.
- Natural Oak laminate base and frame.
- Includes (11) pegboard safety hooks.
- 2" diameter swiveling easy-rolling casters.

MODEL	WIDTH	DEPTH	HEIGHT
5570	24"	18"	57"

Model 5570-100 Accessorized Weight Rack



MODEL	CONSISTS OF
5570-100	(1) 5570 rack, (1) 5581 Set of (16) Hugger® Weights

B. Model 5565 "Combo" Weight/Dumbbell Mobile Rack



FEATURES:

- Almond pegboard panel on front of rack for storage of optional Hugger® Weights, includes (11) hooks.
- Black laminate panel on back equipped with hooks for (20) dumbbells.
- Oak laminate frame & base.
- 3" diameter swiveling, easy-rolling casters.

MODEL	WIDTH	DEPTH	HEIGHT
5565	24"	18"	57"

Model 5565-100 Accessorized Combo Rack



MODEL	CONSISTS OF
5565-100	5565 Rack, (1) 5581 Set of (16) Hugger® Weights, (2) 5505 Sets of (10) Dumbbells each

Model 5571 Pegboard Hooks



- Set of (11) additional Pegboard Hooks with ball tip ends for safety.

MODEL	LENGTH	QUANTITY
5571	4"	11

Common Features

LATEX-FREE RESISTIVE EXERCISE BANDS:



- For multi-functional exercise use.
- Bands are 50 yds. long x 4" wide. Supplied in dispensing box.
- Available in (5) color-coded resistance levels.
- Non-slip, non-latex elastomer.
- No powder, no odor.
- Use this clinically tested allergen-free safer alternative.

C. Model 5531 Set of (5) REP Bands™

Consists of the following:

LEVEL	COLOR
Extra-Light	Peach
Light	Orange
Medium	Green
Heavy	Blue
Extra-Heavy	Plum

REP Band
Resistive Exercise Products

C. Model 5531

