

A. Model 6020-709



B. Model 6040-709



Common Features

TILT TABLES — MODELS 6020-709, 6040-709, 6056-709:



- Gray vinyl upholstery.
- 2" High-Density urethane foam top.
- Tilt angle indicator from 0° to 90°.
- Three extra-wide safety straps.
- Set of (4) 4" locking casters.
- Heavy duty steel frame is powder coated in neutral cream finish.
- 3½" clearance under frame.*
- Designed for rehab & therapy.
- Not designed for cardiac testing.

A. Model 6020-709 Manual Crank Tilt Table



- Manual crank tilt mechanism locks in at any angle. Crank handle is removable.
- Fixed table height is 33" high.
- Removable padded footboard.

MODEL	LENGTH	WIDTH	HEIGHT
6020-709	78"	28"	33"

OPTIONS

FOR MODELS 6020-709, 6040-709, 6056-709

- 6026** – Additional Safety Strap
- *077** – (4) 6" casters increase height and clearance by 2"

B. Model 6040-709 Electric Tilt Table



- Heavy duty electrical tilt mechanism.
- Pneumatic hand control operable from either side of table.
- Fixed table height is 33" high.
- Removable padded footboard.
- Tilts from horizontal to 90° at fixed rate of speed of 3° per second.
- Limit switches.
- Hospital grade plug.
- 110V, 60 Hz, 3 Amps.

MODEL	LENGTH	WIDTH	HEIGHT
6040-709	78"	28"	33"

C. Model 6056-709



C. Model 6056-709 Electric Hi-Lo Tilt Table

SPECIAL FEATURES:

- (2) Heavy duty motors independently operate height and tilt.
- Electrical Hi-Lo motorized lift from 25" to 39" high.
- Electrical tilt mechanism 0° to 90°.
- Pneumatic hand control is accessible from either side of the table.
- (2) Patient support handles.
- Stainless steel fixed footboard.
- Limit switches.
- Hospital grade plug.
- 110V, 60 Hz, 3 Amps.

MODEL	LENGTH	WIDTH	HEIGHT
6056-709	78"	28"	25"- 39"