

Model 6058 elevates to 33" high



Model 6058 lowers to 19" for wheelchair transfer



A. Model 6058

**A. Model 6058
Wheelchair Accessible
Hi-Lo Tilt Table**



FEATURES:

- Unique wheelchair accessible design is ideal for rehabilitation or nursing home use.
- Electric operated actuators provide quick and smooth control of the Tilt and Height adjustments.
- Electric Hi-Lo lift mechanism adjusts from 19" low wheelchair height to 33" high.
- Tilts from horizontal to 90°. Unit must be raised to maximum height before tilting.
- Easy-to-read tilt angle indicator on both sides.
- Pivoting Footboard with rubber tread. Angle of the Footboard tilt is adjustable with spring loaded dual controls.
- "Hands-free" Foot Control operation.
- Includes unique Emergency Tilt Release Handle that manually lowers tilted table back to horizontal position.
- Comfortable 76" x 26" cushioned top is available in Dove Gray Vinyl.
- (3) safety straps with "hook and loop" closures.
- (4) 3" dia. swiveling locking casters.
- Heavy-duty steel frame is cream powder coated epoxy finish.
- 110V, 60Hz. 3.5 amps. Hospital Grade Plug.
- Not designed for cardiac testing use.

MODEL	LENGTH	WIDTH	HEIGHT
6058	76"	26"	19"-33"

**B. Model 6045-709
Bariatric Electric Tilt Table**

FEATURES:



- Electric Tilt Mechanism with heavy-duty motor for up to 500 lb. patients.
- Tilt indicator displays angle of tilt 0° to 90°.
- Extra-wide 34" x 77" long top upholstered in 709 Gray vinyl with 2" thick High-Density urethane foam.
- 23" extra-wide x 18" deep fixed foot board.
- Three extra-long patient safety straps.
- Fixed table height of 33".
- Hand Control operable from either side of the table.
- Set of (4) 3" locking casters.
- Heavy-duty gauge steel frame is powder coated in neutral cream finish.
- 3.5" clearance under frame.
- 110V, 60 Hz, 3.5 Amps.
- Hospital grade plug.

MODEL	LENGTH	WIDTH	HEIGHT
6045-709	77"	34"	33"



B. Model 6045-709