

**A. MODEL 6071**



**B. MODEL 6074**



## COMMON FEATURES

HI-LO SECTIONAL TABLES –  
MODELS 6071, 6074:



- “Hands-free” foot control operation.
- Open base allows closer access to patient.
- Powder coated heavy duty steel frame
- Quiet and stable, single column actuator adjusts from 24” to 33” H.
- Single column actuator saves space, no top shift when elevating.
- Hospital grade power cord.
- 100V-240V, 50/60Hz, 3.15 Amps

### A. MODEL 6071 2-SECTION HI-LO TREATMENT TABLE WITH ARMRESTS AND ADJUSTABLE GLIDES

*SPECIAL FEATURES:*

- 2-Section Power Table can be used for a wide variety of purposes.
- Folding armrests for extra patient comfort.
- Includes paper dispenser (factory installed).
- 2-Section top is 75” L x 28” W with 2” thick high density foam.
- Head/Backrest is 25” L x 28” W & elevates up to 67° with dual release friction lock.
- Adjustable glides allow for easy leveling.
- Agate Blue vinyl upholstery.

### B. MODEL 6074 3-SECTION MOBILE HI-LO TREATMENT TABLE WITH LOCKING CASTERS



*SPECIAL FEATURES:*

- 3-Section Power Table can be used for a wide variety of purposes.
- (4) 3” diameter, swiveling and locking casters.
- 3-Section top is 78” L x 28” W with 2” thick high density foam.
- Contoured face/nose cutout in head section.
- Head section, 19” L x 28” W adjusts from flat to 90° below horizontal.
- Middle section, 21” L x 28” W is stationary.
- Leg/Backrest section, 38” L x 28” W, elevates up to 70° with dual release friction lock.
- Agate Blue vinyl upholstery.

| MODEL | LENGTH | WIDTH | HEIGHT  |
|-------|--------|-------|---------|
| 6071  | 75”    | 28”   | 24”-33” |

| MODEL | LENGTH | WIDTH | HEIGHT  |
|-------|--------|-------|---------|
| 6074  | 78”    | 28”   | 24”-33” |

**C. MODEL 6053**



### C. MODEL 6053 3-SECTION HI-LO MULTI-POSITION TREATMENT TABLE

- Multi-position table can be used for a wide variety of therapy and treatment functions.
- 76”L x 27”W 3-section 1 ½” thick high density gray vinyl top features contoured face/nose cutout in head section.
- Motorized height adjustment from 19” wheelchair height to 37”H.
- “Hands-free” foot control operation.
- Adjustable head section, 15” x 27”, raises up to 45° and lowers down to 90° below horizontal.



- 22” L x 27”W Middle seat section is stationary
- 39”L x 27”W Leg/Backrest section raises easily up to 73° with sturdy friction locking device and dual release handles.
- Top sections have welded steel support frames.
- Includes (4) 3” diameter swiveling casters with locking brakes.
- 110V, 60 Hz, 1.8 amps. U.L. Listed.

| MODEL | LENGTH | WIDTH | HEIGHT  |
|-------|--------|-------|---------|
| 6053  | 76”    | 27”   | 19”-37” |