




HAUSMANN WORK-WELL™ SYSTEM



System III

WORK-WELL™ SYSTEM I: Start-Up Package

- Natural Oak high pressure laminate protective finish.

(7) Components as listed:

- 8901 - Lift and Place Rack Station
- 8911 - Carry-All Weight Box
- 8912 - Tool Kit Weight Box
- 8913 - Stockroom Crate Weight Box
- 8914 - Packing Carton Weight Box
- 8952 - Set of 16 Disc Weights
- 2216 - Footstool

WORK-WELL™ SYSTEM II: Intermediate Package

Same as System I, with the addition of:

- 8930 - Sled Station

WORK-WELL™ SYSTEM III: Extended-Capacity Package

Same as System II, with the addition of:

- 8920 - Shovel Station, Shovel, Sand
- 8940 - Work Table
- 8950 - Ergo Chair

MODEL 8930 PUSH-PULL SLED STATION

- Simulates tasks that involve moving and conveying equipment and materials.
- (4) Removable weight pins for optional disc weights.
- Self-lubricating polyethylene sled runners.
- 24" wide x 40" deep x 46" high. 75 lbs. weight capacity.

MODEL 8901 LIFT AND PLACE RACK STATION


- Simulates lifting and loading tasks.
- (2) Adjustable shelves "lock" into side mounted support rails.
- 30" wide x 18" deep x 78" high. 100 lbs. weight capacity.

MODEL 8920 SHOVEL AND LEVERING STATION

- Simulates tasks that involve shoveling, levering and lifting.
- (3) Fixed "gravity flow" angled shelves. Slip-resistant treads.
- Includes scoop shovel, shovel hook and 100 lbs. of sand.
- 36" wide x 72" deep x 61" high.

MODEL 8952 DISC WEIGHT SET OF 16

- (4) 2 1/2 lb., (8) 5 lb., (4) 10 lb. discs.



Model 8985

MODEL 8985 STRENGTH TEST PLATFORM

- 24" x 24" portable anti-slip wood platform with hook and 36" long chain. Includes carrying handle.



Model 8925

MODEL 8925 MULTI PURPOSE BOLT MODULE

- Bolt board can be used either horizontally, vertically, at 45° angle or in #8901 rack.
- Includes (30) various bolts and nuts in (5) different types and sizes. Includes (2) screwdrivers, (2) wrenches, (1) Allen wrench.



Model 8911

Model 8913

Model 8912

Model 8914

WEIGHT BOXES

Common Features:

- (2) Hand grip slots on each side. Durable dowel construction.
- (1) Removable weight pin for optional disc weights.
- 3/4" plywood covered with Oak high pressure laminate.

MODEL 8911 "CARRY-ALL" WEIGHT BOX

- 11 1/2 lbs. empty. 14" wide x 14" deep x 14" high.

MODEL 8912 "TOOL KIT" WEIGHT BOX

- 16 lbs. empty. 19" wide x 14" deep x 14" high.

MODEL 8913 "STOCKROOM CRATE" WEIGHT BOX

- 16 lbs. empty. 17" wide x 17" deep x 14" high.

MODEL 8914 "PACKING CARTON" WEIGHT BOX

- 17 lbs. empty. Hinged lid. 14" wide x 14" deep x 14" high.