

Common Features

SMART CARTS™ & BEDSIDE CABINETS—
MODELS 9018, 9029, 9038:



- Top laminate surface: 17" W x 17½" D.
- All cabinets have 3mm high impact edging.
- Easy clean polystyrene Drawer interiors with rounded inside corners.
- 4" Satin Finish wire drawer pulls with Oak cabinet, or Black nylon wire pulls with Gray cabinet.
- 2½" twin wheel hooded swivel Nylon Casters on all Mobile models.
- Ships fully assembled.

COLORS:

Laminate – Natural Oak or Folkstone Gray.
Optional 5904 Wild Cherry & 7909 Fusion Maple (9038-20 Only)

A. Model 9038-20-346



Close up of 3mm edging.

Model 9018-20 Mobile Cabinet (not shown)

- (1) 5½" drawer & (1) hinged door.
- Includes (1) 11½" deep adjustable shelf.

MODEL	LENGTH	WIDTH	HEIGHT
9018-20-346	18½"	18½"	30"
9018-20-927	18½"	18½"	30"

B. Model 9018-10 Stationary Bedside Cabinet

- (1) 5½" Drawer & (1) hinged door.
- Use as bedside or storage cabinet.
- Includes (1) 11½" deep adjustable shelf.

MODEL	LENGTH	WIDTH	HEIGHT
9018-10-346	18½"	18½"	28"
9018-10-927	18½"	18½"	28"

B. Model 9018-10-346



A. Model 9038-20 Mobile Cabinet

- (2) 5½" Drawers & (1) 8½" deep Drawer.

MODEL	LENGTH	WIDTH	HEIGHT
9038-20-346	18½"	18½"	30"
9038-20-927	18½"	18½"	30"
9038-20-5904	18½"	18½"	30"
9038-20-7909	18½"	18½"	30"

C. Model 9029-24 Mobile Cart

- Accommodates EKG or compact modality equipment.
- Available only in Folkstone Gray.
- Drop lid Door & open storage.

MODEL	LENGTH	WIDTH	HEIGHT
9029-24-927	18½"	18½"	30"

C. Model 9029-24-927



**CARTS & CABINETS
NOW STANDARD WITH
3MM EDGING**

- High impact edging for increased product durability.
- All edges are trimmed and buffed which eliminates sharp corners.
- All finished edges are 3mm and will not peel.




Scan for more info on 3mm edging

D. Model 9029-8 Mobile Cabinet

- (1) 5½" Drawer
- Open storage.

MODEL	LENGTH	WIDTH	HEIGHT
9029-8-346	18½"	18½"	30"
9029-8-927	18½"	18½"	30"

D. Model 9029-8-927



See page 90 for additional wider/taller carts